

ACTINIA[®] Hip Stem

MDR Certification



Primary Hip Stem System

ACTINIA® includes a complete range of cementless and cemented hip stem implants. The double-tapered stem shape ensures good primary stability. The flattened neck geometry increases the range of motion.



Proven Design

- on the market since 2012
- more than 60,000 implantations



Size Spectrum

- 11 sizes standard
- 10 sizes lateralised



Broach Handles


- straight
- GIS® (incl. WP-version)
- double offset
- DAA GIS®



Material

- cless.: implatan® TiAl₆V₄
- cem.: implavit® CoCrMo





walk with us

implantcast GmbH

Lüneburger Schanze 26
21614 Buxtehude
Germany

phone: +49 4161 744-0

fax: +49 4161 744-200

info@implantcast.de

www.implantcast.de

CE 0482



ACTNIE-220916

