

BICANA[®] Hip Stem



Primary / Revision Hip Stem System

BICANA® includes a complete range of cemented hip stem implants for primary or revision cases. The extensive collar allows physiological proximal load transmission and has a resection angle of 45° to the femoral axis.



Proven Design

- on the market since 2008
- more than 27,700 implantations



Broach Handles

- straight
- double offset
- GIS®




Size Spectrum

- Primary 8 resp. 9 sizes
- Revision 126°: 9 sizes / 135°: 6 sizes



Material

- cem.: implavit® CoCrMo



walk with us

implantcast GmbH

Lüneburger Schanze 26
21614 Buxtehude
Germany

phone: +49 4161 744-0

fax: +49 4161 744-200

info@implantcast.de

www.implantcast.de

CE 0482

