

ECOFIT[®] Hip Stem



Primary Hip Stem System

ECOFIT® includes a complete range of cementless and cemented hip stem implants. Cementless short stems are also part of the portfolio. Different CCD angles and coating variants offer a variety of treatment options.



Proven design

- on the market since 2004
- more than 220,000 implantations



Broach handles

- straight
- **GIS®** (incl. WP version)
- double offset
- DAA **GIS®**




Short stem


- version of the **ECOFIT®** hip stem 133° shortened by 35 mm



Material

- cless.: implatan® TiAl₆V₄
- cem.: implavit® CoCrMo





Surface
cem.: highly polished
matt blasted

Coating
class.: implaFix® cpTi
implaFix® Duo, cpTi + HA
cem.: TiN (optional for 138° stem)

CCD angle
138° std. + lat.
133° std. + lat.
123° std.
lateralisation without change
of the CCD angle

ic taper
12/14 mm

walk with us

implantcast GmbH

Lüneburger Schanze 26
21614 Buxtehude
Germany

tel.: +49 4161 744-0
fax: +49 4161 744-200

info@implantcast.de
www.implantcast.de

CE 0482

