


Full Pelvis Solution



MUTARS® PRS



indication

acetabular bone defects



anchorage

equatorial PressFit and
boneside **EPORE®** structure



screw options

9 bohre holes for spongiosa
screw flat head Ø 6.5mm
and 3 bore holes for
spongiosa screw Ø 8mm



size range

52 - 72 mm in 4 mm steps



MUTARS®

PRS reinforcement cage



indikation

pronounced acetabular
bone defects



compatibility

with **MUTARS®** PRS




variants

short (5 bore holes) and long
(10 bore holes) proximal flap



screw options

bore holes for spongiosa
screw flat head Ø 6.5mm



EPORE[®] buttress & shim augments



indication

acetabular bone defects



compatibility

with all ic-hip cups

(only **EPORE[®]** buttress augments)



screw options

bore holes for spongiosa

screw flat head Ø 6.5 mm



size range

EPORE[®] buttress augment:

4 sizes; **EPORE[®]** shim augment:

3 sizes



ic-restrictor acetabulum



indication

acetabular bone defects



design feature

thread for a guided setting
of the implant




compatibility

with all ic-hip cups



size range

3 sizes (26, 32, 38 mm)



Full Pelvis Solution


EPORE[®]

acetabulum spacer

8 mm



removable
crunch plugs



indication

acetabular bone defects



compatibility

with all ic-hip cups





screw options

bore holes for spongiosa screw
flat heads Ø 6.5mm and
cancellous screw Ø 8mm



size range

26 sizes



aim

conversion of an uncontained defect into a contained defect



modular acetabular reamer guide

for preparation of **EPORE®** acetabular spacer



combination

possibility for adjustment of the acetabular spacer heights (10, 15, 20)



size range

modular acetabular reamer guides are available in sizes 52, 56, 60, 64 and 68mm

walk with us

implantcast GmbH

Lüneburger Schanze 26
21614 Buxtehude
Germany

tel.: +49 4161 744-0
fax: +49 4161 744-200

info@implantcast.de
www.implantcast.de

CE 0482



Document number:
FUPEFLEN-201223